

Respirator Fit Testing

OSHA Compliance & Safety

What % of workers fail to achieve adequate facial seal of N95 masks?

50%



Respirator fit testing ensures that employees are safe from breathing in harmful contaminants and pathogens from their environment.



Online Booking



Results provided to patient



Aggregate and Individual Data for Employers



Efficient processing of large groups

Program Overview

We provide qualitative respirator fit testing to help employers comply with OSHA requirements. Qualitative fit testing is a pass/fail test method that uses your sense of taste or smell, or your reaction to an irritant in order to detect leakage into the respirator facepiece.

Quantitative fit testing uses a machine to measure the actual amount of leakage into the facepiece. We can also provide this service upon request through our occupational medicine partners.

Program Details

Length: 5 minutes

Space Requirements: 6X6

Min Appts: 40

Coverage: Direct Pay

Book Appt

Provide Medical History

Sign Consent Forms

Have Onsite Experience

Receive Data